

VetoRapid


FOR THE EASY
IDENTIFICATION
OF BACTERIA



VetoRapid

Diagnostic tool to aid in the rapid
detection of mastitis in lactating cows
Narzędzie diagnostyczne do szybkiego
wykrywania mastitis u krów mlecznych

05 | 
TESTS



03 A diagnostic tool to identify mastitis causing agents in dairy cows.

- **Contents:** 5 VetoRapid tests
- **Storage:** refrigerate (2-8°C) and keep away from light. Keep out of the sight and reach of children.

How to use:

- 1 Collect a sterile sample of milk from the affected quarter and label the sample.
- 2 The sample must be kept in the refrigerator at a maximum temperature of 7°C for no longer than 24 hours. Agar must stay dry in the possible presence of contamination plus plate within the incubator until dry.
- 3 Prepare your VetoRapid plate in a clean place and use a sterile needle to guide the milk sample.
- 4 Seal the seal from the milk sample in a cap cap pattern on each part of the plate.
- 5 Load your VetoRapid plate and place it upside down flat at the bottom within the incubator at 25-37°C.
- 6 Only one sample per VetoRapid plate.
- 7 Leave to incubate for 24 hours.
- 8 Identify the bacteria. See QR code on lid.

Danabac in UK: Vetoquinol UK Limited,
Plym Hill Business Park, Dorchester,
Northwest, Wiltshire, UK

08/PL/00016 VetoRapid

ACHIEVE
BETTER MILK
TOGETHER




VetoRapid

Vetorapid is a diagnostic tool that has been designed to aid in the easy identification of the most common bacteria found in mastitic milk samples. It is not intended to replace a high quality diagnostic laboratory but is a useful first line tool in basic identification of the main pathogens that cause clinical mastitis; *E.coli*, Staphylococci and Streptococci¹.

A university of Glasgow trial found that Vetorapid provided a rapid preliminary identification of five common causes of bovine mastitis (*Escherichia coli*, *Staphylococcus aureus*, coagulase-negative staphylococci, *Streptococcus uberis* and *Enterococcus* species) under field conditions².


Introducing Vetorapid

Vetorapid dish compartments for bacterial identification:




Section 1

Selective for Gram negative bacteria



Section 2

Selective for staphylococci



Section 3

Selective for streptococci and enterococci

How to correctly take a sample for Vetorapid testing

1. Take a clean milk sample from the affected quarter - See milk sampling procedure below for details
2. Place the sample in the fridge for storage
3. Prepare your plate in a clean area and take a sterile swab or a loop to dip into the milk sample
4. Spread the milk sample on the agar plate, using the zig-zag technique
5. Label the plate and incubate at 35-37°C with the lid on (lid must be facing downwards)
6. Do the next sample the same way, with a new plate. Only one sample should be analysed per plate
7. Results within 18-36 hours*
8. Follow the diagnostic table, sector by sector.

NB. Plates should be dry. If there is condensation on the agar, let it dry in the incubator before use.

Taking a milk sample on farm:



How to apply the milk sample to the Vetorapid plate

We recommend the use of a loop (metal loop or disposable loop) to apply the milk sample to the Vetorapid plate as results are easier to analyse with this method.

- single colonies are easier to detect and there is better differentiation between possible growth and cloudiness caused by the smeared milk
- plate the milk as shown in the picture **A**. Do this for all three sections.


[N.B. Sterile cotton swabs can also be used but can cause cloudiness making it harder to differentiate between this and possible growth; however this method can improve sensitivity].

Incubation of Vetorapid

Once the sample has been applied to the Vetorapid plate, label and incubate it at 35-37°C with the lid on facing downwards.

Results should show within 18-36 hours*.
Most colonies are detectable within 24 hours.

*Some cultures can take up to 48 hours.




Veto**Rapid**

How to read the results by plate sector

Gram negative coliform agar (Sector 1)


Generally *E. coli* and other coliforms are clearly detectable in less than 18 hours.



***E. Coli*;**

Colony colour: Dark blue-violet

Agar colour: No change in agar colour



Other coliforms,

example: Klebsiella;


Colony colour: Red-purple

Agar colour: No change in agar colour

How to read the results by plate sector

Staphylococcus agar (Sector 2)


Generally staphylococci (e.g. *Staph. aureus* and coagulase negative staphylococci) are clearly detectable after 24 hours.



***Staph. aureus*;**

Colony colour: Yellow-golden

Agar colour: Gold



Coagulase-Negative Staph. (CNS),

example: Staph. epidermidis;


Colony colour: Clear

Agar colour: No change in agar colour

How to read the results by plate sector

Streptococcus and enterococcus agar (sector 3) - agar aesculin positive

Aesculin positive bacteria are clearly detectable by the black discolouration of the agar after 24 h; of the aesculin positive bacteria, *Strep. uberis* accounts for $\geq 95\%$ of cases and enterococci $\leq 5\%$ of cases.




Aesculin positive streptococci,

example: *Strep. uberis*;

Colony colour: Black

Agar colour: Black



Aesculin positive enterococci;

Colony colour: Black
(Yellow in
sector 2)


Agar colour: Black

How to differentiate between *Strep. uberis* and enterococci:


- Generally enterococci and *Strep. uberis* cannot be differentiated by this mastitis test (a salt tolerance test is necessary).
- Some enterococci strains will also grow slightly on the Sector 2 agar. Look for black discolouration of the Sector 3 agar after 24 hours, combined with a yellow discolouration of the Sector 2 agar after 48 hours (less obvious than the yellow discolouration caused by *Staph. aureus*).

How to read the results by plate sector

Streptococcus and enterococcus agar (sector 3) - agar aesculin negative



In Sector 3 agar it is difficult to differentiate between α - and β - haemolysis after using a swab, we recommend using loops. Smaller colonies of aesculin negative streptococci are difficult to detect before 24 hours but are usually more visible after 48 hours.



Aesculin negative *Strep. dysgalactiae*;

Colony colour: Clear

Agar colour: Red-brown

α -haemolysis: Green coloured corona around the colonies.

Aesculin negative *Strep. agalactiae*;






Colony colour: Clear

Agar colour: Red-brown

β -haemolysis: Clear bright corona around the colonies

Summary of all streptococcus and enterococcus agar readings:





VetoRapid

*Some cultures can take up to 48 hours.

ACHIEVE BETTER MILK TOGETHER



VetoRapid

References: 1. Bradley A.J, Leach K.A, Breen J.E, Green L.E, Green M.J. Survey of the incidence and aetiology of mastitis on dairy farms in England and Wales. 2007 Veterinary Record 160, 253-258. 2. Viora L, Graham E.M, Mellor D.J, Reynolds K, Simoes P.B.A, Geraghty, T.E. Evaluation of a culture-based pathogen identification kit for bacterial causes of bovine mastitis. 2014 Veterinary record 175, 89.

Further information is available on request from: Vetoquinol UK Limited, Vetoquinol House, Great Slade, Buckingham Industrial Park, Buckingham, MK18 1PA. UK: Tel: 01280 814500 Fax: 01280 825460 ROI: Tel: 1800 406117 Fax: 1800 406116 Email: office@vetoquinol.co.uk Website: www.vetoquinol.co.uk

Please use medicines responsibly. For further information please visit www.noah.co.uk/responsible

